

**Product Name :**  
Pathology Lab Semi Auto Analyzer**Product Code :**  
MED-PTH0016**Description :**

Pathology Lab Semi Auto Analyzer

**Technical Specification :**

The Semi-Auto Analyzer helps in measuring chemical indexes by examining blood and other body fluid. Further, it is combined with other clinical information, to assist in evaluating organs function, diagnosing diseases, recognizing disease gene and determining the norm for future therapy. The in-built thermal printer of the analyzer prints the results and makes, diagnosing a disease simpler process.

**Product Details :**

Application: Laboratory Use  
Display Type: TFT touch-screen  
Display Size: 7.0 inch  
Number of Test Modes: 2  
Test Mode Type: Flowcell or cuvette  
Packaging Type: Boxes

**SALIENT FEATURES**

7.0" TFT touch-screen & pop-up keypad  
Two test mode: flowcell or cuvette  
Supports bi-chromatic tests for end point, fixed-time, kinetics methods  
External keyboard or mouse support via USB  
User-friendly interface and multi-language operation software  
Built-in thermal printer

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,