

**Product Name :**  
Pathology Lab Elisa Washer**Product Code :**  
MED-PTH0017**Description :**

Pathology Lab Elisa Washer

**Technical Specification :**

Elisa Washer is perfect for providing gentle washing of cell based assays with adherent or weakly adherent cells and ensures excellent operational life and easy maintenance. Fully automatic, the washer has easy window operation system and thereby ensures simple functioning. The range is absolutely compatible with U or V, flat bottom plates and is has variable soaking and shaking time.

**Product Details :**

Application: Clinical

Wash Modes: 2

Wash Mode Type: Interchangeable manifolds for 8- and 12-well strips

Bottom Detection: Automatic

Compatibility: With flat, U- or V-bottom plates

Packaging Type: Boxes

**FEATURES**

Two rinsing bottles avoids bottle change

Automatic bottom detection and collision protection

Variable shaking and soaking time

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)