

Product Name :
Pathology Lab Elisa Reader**Product Code :**
MED-PTH0018**Description :**

Pathology Lab Elisa Reader

Technical Specification :

Elisa Reader offered by us is extensively used for agglutination reading. This reader is provided with data analysis software and ensures simple operations. The reader has 7.8" touch-screen with pop-up keypad and touch pen, which allows easy readings. In addition to this, the reader thermostatically balanced in nature and made available at convenient prices.

Product Details :

Application: Hospital

Display Type: LCD

Packaging Type: Boxes

Display Size: 7.8 in touch-screen

Memory Storage: Up to 20,000 results

Interface: RS232

FEATURES

In-Built thermal printer

Variable speed shaking and time modes

Superior -capacity memory stores up to 20,000 results

7.8" touch-screen with touch pen and pop-up keypad

Influential data management software for computer installation (optional)

RS232 for external PC control, 2 USB ports for keyboard, external mouse, and printer

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,