


Product Name :
Pathology Lab Blood Cell Counter**Product Code :**
MED-PTH0019**Description :**

Pathology Lab Blood Cell Counter

Technical Specification :

Blood Cell Counter is completely automatic hematology analyzer with 19 parameters for microsampling technology and CBC test. Compact in size, the counter is generally used for counting the RBC in the blood. Further, the user-friendly interface of the counter allows simple and effective workflow. This type of counter ensures accurate results and efficient performance.

Product Details :

Application: Hospital, Clinical
Operation Mode: Semi Automatic
Sample Volume: 13 ul (whole blood)
Storage Up to: 20,000 sample results
Throughput: 30 samples per hour
Number of Counting Modes: 2
Counting Mode Type: Whole blood and Prediluted

FEATURES

Three-part differentiation of WBC, 19 parameters +3 histograms
Large Color LCD Display

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,