


**Product Name :**  
Pathology Lab Blood Cell Counter**Product Code :**  
MED-PTH0019**Description :**

Pathology Lab Blood Cell Counter

**Technical Specification :**

Blood Cell Counter is completely automatic hematology analyzer with 19 parameters for microsampling technology and CBC test. Compact in size, the counter is generally used for counting the RBC in the blood. Further, the user-friendly interface of the counter allows simple and effective workflow. This type of counter ensures accurate results and efficient performance.

**Product Details :**

Application: Hospital, Clinical  
Operation Mode: Semi Automatic  
Sample Volume: 13 ul (whole blood)  
Storage Up to: 20,000 sample results  
Throughput: 30 samples per hour  
Number of Counting Modes: 2  
Counting Mode Type: Whole blood and Prediluted

**FEATURES**

Three-part differentiation of WBC, 19 parameters +3 histograms  
Large Color LCD Display

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)