

Product Name :
Pathology Lab Microscope**Product Code :**
MED-PTH0020**Description :**

Pathology Lab Microscope

Technical Specification :

Microscope offered by us has quadruple revolving nosepiece and allows interchangeability of optics. Helps to clearly see the micro organisms, this is provided with mechanical stage right hand drive. The in-built 6v, 200W halogen light illuminator, in the microscope ensures high performance.

Product Details :

Application: Clinical, Laboratory
Surface Illumination: Halogen
Halogen Voltage Required: 6 V
Halogen Wattage: 200 W
Stage Length: 120 mm
Stage Width: 132 mm

FEATURES

Antifungus optics (sealed optical parts)
Interchangeability of optics
Abbe condenser & light relay system are fitted with high performance aspheric lenses for bright & crisp image
Illumination system through SMPS circuit etc.
Coaxial coarse and fine roller guide (Comfortable position of focusing knobs). Tension adjustment on right
Quadruple revolving nosepiece
With mechanical stage right hand drive

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,