


Product Name :
Pathology Lab Centrifuge**Product Code :**
MED-PTH0021**Description :**

Pathology Lab Centrifuge

Technical Specification :

The Centrifuge operates under an irregular filtration process and is perfect to the competent separation of pharmaceutical products and bulk chemicals. This type of centrifuge is used along with tiny batch sizes and where products are changed often with less process costs. The digital speed indicator and Stepless speed regulator, makes this product highly efficient.

Product Details :

Application: Hospital, Clinical
Basket Material: SS
Display Type: Digital
Frequency: 50 / 60 Hz
Power Source: Battery
Packaging Type: Boxes

SALIENT FEATURES

Digital Pace Indicator
Stepless speed controller
0-60 minutes digital countdown timer
Press Switch for quick acceleration to full speed

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,