


**Product Name :**  
Pathology Lab Centrifuge**Product Code :**  
MED-PTH0021**Description :**

Pathology Lab Centrifuge

**Technical Specification :**

The Centrifuge operates under an irregular filtration process and is perfect to the competent separation of pharmaceutical products and bulk chemicals. This type of centrifuge is used along with tiny batch sizes and where products are changed often with less process costs. The digital speed indicator and Stepless speed regulator, makes this product highly efficient.

**Product Details :**

Application: Hospital, Clinical  
Basket Material: SS  
Display Type: Digital  
Frequency: 50 / 60 Hz  
Power Source: Battery  
Packaging Type: Boxes

**SALIENT FEATURES**

Digital Pace Indicator  
Stepless speed controller  
0-60 minutes digital countdown timer  
Press Switch for quick acceleration to full speed

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)