

Product Name :
Pathology Lab Bpm Blood Pressure Monitor**Product Code :**
MED-PTH0023**Description :**

Pathology Lab Bpm Blood Pressure Monitor

Technical Specification :

The BPM automatic blood pressure monitor features a motion detector if it is moving during the pressure measurement (Movement Detector). It works with Real Fuzzy Technology. Hypertension Risk Detector Display the average of 3 recent measurements, save up to 90 data, easy to use, just one press.

Product Details :

Pressure: 30 - 260 mm

Mercury Precision: +/- 3 mm

Mercury Pulse: 40 - 199 beats / min +/- 5% accuracy

Memory: 90 values

Adaptor: Power 6 V