

**Product Name :**  
Pathology Lab Pulse Oximeter

**Product Code :**  
MED-PTH0028



**Description :**

Pathology Lab Pulse Oximeter

**Technical Specification :**

Display Mode: 3 DISPLAY MODES  
Oxygen Saturation Range: 25-100%  
Parameters: Pulse Rate, PULSE BAR  
Pulse Rate Measuremet Accuracy: 30 - 99BPM  
Pulse Rate Measuremet Range: 30 - 235BPM  
SPO2 Measuremet Range: 60% - 100%