

Product Name :
Pathology Lab Urine Analyzer**Product Code :**
MED-PTH0034**Description :**

Pathology Lab Urine Analyzer

Technical Specification :

LCD display, keyboard operation.

Save data automatically.

Date and number can be input.

The colour and limpidness can be input.

Single test and continuous test are available.

Adopt luminescence source with high brightness to minimize interference of the environment light.

Throughput: 120Tests/hour

Display: LCDdisplay

Test Items: PH, Nitrite, Protein,Glucose, Occult blood, Bilirubin,Urobilinogen, Ketone, Specific gravity, Leukocytes,Ascorbic Acid,Color,Limpidness

Principle: Reflection of high-brightness cold light

Printer: Built-in thermal printer,paper type 57mm

Storage: 1000 patient datas

Working Environment: Temperature:0~40?; humidity: 10~80%

Interface: Standard RS-232port

Power Supply: AC220±10%, 50/60Hz, 110±10%,60Hz

Medical Hospital Equipments Manufacturers

Website: [Medical Hospital Equipments Manufacturers](#),

NAUGRA[®]