

**Product Name :**  
Pathology Lab HB1Ac Machine

**Product Code :**  
MED-PTH0036



**Description :**

Pathology Lab HB1Ac Machine

**Technical Specification :**

Sample type: Venous blood, Peripheral finger blood  
Sample volume: 5?l  
Sampling mode: Full-auto 25-position sample turntable  
Testing mode: Auto test for batch  
Testing item: Glycated Hemoglobin HbA1c  
Display: 320x240 LCD Graphic Display  
Data storage: 1000 test report  
Power supply: AC220V±22V 50Hz 65VA