

**Product Name :**  
Pathology Lab Eliza Plate Washer and Reader

**Product Code :**  
MED-PTH0037



**Description :**

Pathology Lab Eliza Plate Washer and Reader

**Technical Specification :**

8-channel vertical optical system.  
Super fast reading speed of 7 seconds.  
10.1-inch Large Touch Screen with visual layout function, what you see is what you get.  
Hardware: ARM architecture high-performance low-power processor, with ADI high-precision AD conversion chip, the measurement is more accurate and stable.  
Software: Self-developed core algorithm, with independent intellectual property rights, flexible use of multiple calculation methods.