

Product Name :
Pathology Lab Eliza Plate Washer and Reader

Product Code :
MED-PTH0037



Description :

Pathology Lab Eliza Plate Washer and Reader

Technical Specification :

8-channel vertical optical system.
Super fast reading speed of 7 seconds.
10.1-inch Large Touch Screen with visual layout function, what you see is what you get.
Hardware: ARM architecture high-performance low-power processor, with ADI high-precision AD conversion chip, the measurement is more accurate and stable.
Software: Self-developed core algorithm, with independent intellectual property rights, flexible use of multiple calculation methods.