

Product Name :
Pathology Lab Immunoassay Analyzer**Product Code :**
MED-PTH0038**Description :**

Pathology Lab Immunoassay Analyzer

Technical Specification :

Modes of Operation: Random, Batch and STAT
Separation method: Magnetic separation technology
Sample Positions: 60 (Each position could be used as an emergency position)
Reagent Positions: 30 pcs (2-8°C refrigeration)
Incubation Positions: 50
Time to 1st Result: 18 Minutes
Sample Volume: 15-100ul
Detection Precision: CV?5%
Carry-over Rate: ?10-5
Correlation Coefficient: r?0.99
Calibration Stability: 28 Days
Calibration Type: 6/7-point calibration
Software system: Windows 7/8/10, 32 or 64 bit
Bi-direction, support: HL7 protocol
Intuitive User Interface, Intelligent Data Management
Interface: TCP/IP Network interface
Other Function: Liquid level detection
Anti-collision function
Sample barcode scanning
Reagent radio frequency identification
Intelligent alarm prompt

Power Supply: AC110/220V,50/60Hz

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,