

**Product Name :**  
Pathology Lab Paraffin Wax Bath

**Product Code :**  
MED-PTH0043



**Description :**

Pathology Lab Paraffin Wax Bath

**Technical Specification :**

Capacity: 10L  
Temp. Range: RT-100°C adjustable  
Temp. Precision:  $\pm 1^\circ\text{C}$   
Power Supply: AC 220V $\pm 10\%$  50Hz (standard model);  
AC110V $\pm 10\%$ : 60Hz  
Consumption: 1200W