

Product Name :
Pathology Lab Grossing Workstation**Product Code :**
MED-PTH0044**Description :**

Pathology Lab Grossing Workstation

Technical Specification :

External Size (W*D*H): 1800*750*1950mm

Working Area Size (W*D*H): 1695*430*900mm

Work Surface Height: 900mm

Consumption: 300W

Optional water heater: 1500W

Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

Exhaust Airflow Volume: 500-600m³/h

Material Main Body: 304 stainless steel

Work Table: 304 stainless steel

Standard Accessories: Spotlight;

UV lamp;

Fluorescent lamp;

Exhaust system;

4m exhaust duct;

Pipe strap;

Auto flushing device;

Stainless steel hanger;

Adjustable hot and cold water tap;

180mm deep stainless steel water sink, funnel with filter.

Active carbon filter

Optional Accessories: Pulverizer, Electrical water heater, Automatic induction hand dryer, Eye washer, Camera.

Flex magnifier light, laptop bracket

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,