

**Product Name :**  
Pathology Lab Tissue processor**Product Code :**  
MED-PTH0046**Description :**

Pathology Lab Tissue processor

**Technical Specification :**

Fully intelligent design, ensure timely automatic determination and recovery from an abnormal event.  
Manual adjustment can be conducted anytime during the programmed automatic operation; afterwards, the system automatically enters the programmed operation.  
Green inner-cycling air purification system to efficiently absorb and remove poisonous gas.  
Automated fan and light control.  
Internal dry heating mechanism with high-precision temperature control  
Cups Number: 12  
Single Cup Capacity: 1L  
Processing Time Within 0~99 hours for the first cup(BK-TS3D Any length);0-24 hours for othercup.  
Paraffin Cup Temp. Range: RT~99°C  
Temp. Precision: 1°C  
DrippingTime: 0~30 s 0~60s  
Display: LCD  
Agitation Frequency: 2 times/min  
Environment Temp. Range: 0~40°C  
Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**