

Product Name :
Pathology Lab Tissue processor**Product Code :**
MED-PTH0046**Description :**

Pathology Lab Tissue processor

Technical Specification :

Fully intelligent design, ensure timely automatic determination and recovery from an abnormal event.
Manual adjustment can be conducted anytime during the programmed automatic operation; afterwards, the system automatically enters the programmed operation.
Green inner-cycling air purification system to efficiently absorb and remove poisonous gas.
Automated fan and light control.
Internal dry heating mechanism with high-precision temperature control
Cups Number: 12
Single Cup Capacity: 1L
Processing Time Within 0~99 hours for the first cup(BK-TS3D Any length);0-24 hours for othercup.
Paraffin Cup Temp. Range: RT~99°C
Temp. Precision: 1°C
DrippingTime: 0~30 s 0~60s
Display: LCD
Agitation Frequency: 2 times/min
Environment Temp. Range: 0~40°C
Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]