

Product Name :
Pathology Lab Water Bath**Product Code :**
MED-PTH0048**Description :**

Pathology Lab Water Bath

Technical Specification :

Stainless steel inner chamber and upper cover.
Adopt new design structure, which could avoid the leakage of steam.
PID Microprocessor controller, LED display.
Capacity: 3L
Temp. Range: RT.-100?
Temp. Fluctuation: ± 0.5 ?
Temp. Accuracy: ± 1 ?
Time Range: 0~999min
Power: 500W
Voltage: AC110/220VV $\pm 10\%$, 50/60Hz