

Product Name :
Pathology Lab Hot Air Oven

Product Code :
MED-PTH0049



Description :

Pathology Lab Hot Air Oven

Technical Specification :

Capacity: 30L
Temp. Range: RT+10~300°C
Temp. Precision: 0.1°C
Temp. Fluctuation: ±0.5°C
Ambient Temp: 5~40°C
Timing Range: 1~9999min
Shelves No: 2 pcs
Display: LED
Consumption: 800W
Power Supply: AC220V±10%, 50/60Hz;