

Product Name :
Pathology Lab BOD Incubator

Product Code :
MED-PTH0051



Description :

Pathology Lab BOD Incubator

Technical Specification :

LCD touch screen.
Double door design.
Polished stainless steel chamber.
Round angle structure , easy to clean.
Universal castors for easy movement.
CFC-free environmental friendly refrigerant.
Microprocessor PID temperature control: accurate and reliable.
BOD socket integrated inside for additional small equipment.

Safety function:

Temperature protection threshold value can be set.
Auto-off protection and alarm for abnormal situation.
Parameters can be saved automatically when power off or system error.

Website: [SITE](#)

NAUGRA[®]

Medical Hospital Equipments Manufacturers