

**Product Name :**  
Pathology Lab BOD Incubator

**Product Code :**  
MED-PTH0051



**Description :**

Pathology Lab BOD Incubator

**Technical Specification :**

LCD touch screen.  
Double door design.  
Polished stainless steel chamber.  
Round angle structure , easy to clean.  
Universal castors for easy movement.  
CFC-free environmental friendly refrigerant.  
Microprocessor PID temperature control: accurate and reliable.  
BOD socket integrated inside for additional small equipment.

**Safety function:**

Temperature protection threshold value can be set.  
Auto-off protection and alarm for abnormal situation.  
Parameters can be saved automatically when power off or system error.

---

Website: [SITE](#)

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**