

Product Name :
Pathology Lab Humidity Cabinet**Product Code :**
MED-PTH0052**Description :**

Pathology Lab Humidity Cabinet

Technical Specification :

External Size(W*D*H): 700*650*1230 mm
Internal Size(W*D*H): 600*500*540 mm
Tested Opening: Safety Height ? 200 mm(8")
Max Opening: 370mm(15")
Inflow Velocity: 0.53±0.025 m/s
Down Flow Velocity: 0.33±0.025 m/s
HEPA Filter: Two, 99.999% efficient at 0.3?m. Filter life indicator.
Front Window: Motorized. Toughened glass. Anti UV.
Noise: EN12469 ? 58 dB / ? 61 dB
UV Lamp 15W*1: UV timer, UV life indicator, emission of 253.7 nanometers for most efficient decontamination.
LED Lamp: 8W*2
Illumination: ?1000Lux
Consumption: 500W
Socket: Two, total load of two sockets: 500W
Display: LCD display
Control System: Microprocessor
Airflow System: 70% air recirculation, 30% air exhaust
Visual and Audio alarm: Filter replacement, window over height, abnormal air flow velocity.
Material: Work Zone: 304 stainless steel
Main Body: Cold-rolled steel with anti-bacteria powder coating.

Work Surface Height: 750mm with optional base stand
Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#),