

Product Name :
Pathology Lab Humidity Cabinet**Product Code :**
MED-PTH0052**Description :**

Pathology Lab Humidity Cabinet

Technical Specification :

External Size(W*D*H): 700*650*1230 mm

Internal Size(W*D*H): 600*500*540 mm

Tested Opening: Safety Height ? 200 mm(8")

Max Opening: 370mm(15")

Inflow Velocity: 0.53±0.025 m/s

Down Flow Velocity: 0.33±0.025 m/s

HEPA Filter: Two, 99.999% efficient at 0.3?m. Filter life indicator.

Front Window: Motorized. Toughened glass. Anti UV.

Noise: EN12469 ? 58 dB / ? 61 dB

UV Lamp 15W*1: UV timer, UV life indicator, emission of 253.7 nanometers for most efficient decontamination.

LED Lamp: 8W*2

Illumination: ?1000Lux

Consumption: 500W

Socket: Two, total load of two sockets: 500W

Display: LCD display

Control System: Microprocessor

Airflow System: 70% air recirculation, 30% air exhaust

Visual and Audio alarm: Filter replacement, window over height, abnormal air flow velocity.

Material: Work Zone: 304 stainless steel

Main Body: Cold-rolled steel with anti-bacteria powder coating.

Work Surface Height: 750mm with optional base stand
Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#),