

Product Name :
Pathology Lab Laminar Airflow**Product Code :**
MED-PTH0054**Description :**

Pathology Lab Laminar Airflow

Technical Specification :

Air speed adjustable.
Tabletop type,easy to carry and save space.
Microprocessor control system,LED Display.
HEPA Filter with efficiency:99.999% at 0.3 μm.
External Size (W*D*H)mm: 550*460*700
Internal Size (W*D*H): 480*340*370
Display: LED Display
Airflow Velocity: 0.3~0.5m/s
Material: Main Body Cold-rolled steel with anti-bacteria power coating
Work Table: 304 Stainless steel
Pre-Filter: Polyester fiber, washable(for BBS-V500,BBS-V600 and BBS-V700)
Filter: HEPA Filter
Noise: 70dB
Front Window: Manual, 5mm toughened glass, anti UV
Max Opening: 310mm
LED Lamp: 4W*1
UV Lamp: 8W*1 Emission of 253.7 nanometers
Consumption: 100W
Power Supply: AC220V±10%, 50/60Hz; 110V±10%, 60Hz

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,