

**Product Name :**  
Pathology Lab Freezer Refrigerator**Product Code :**  
MED-PTH0055**Description :**

Pathology Lab Freezer Refrigerator

**Technical Specification :**

External Size(W\*D\*H) mm: 500\*550\*1300

Capacity: 130L

Temp. Range: 2°C~8°C

Temp. Accuracy: 0.1°C

Control System: Microprocessor Control, Large Screen LED Display

Alarm: Audible and visual alarm for: High and low temperature, Power failure alarm,

Sensor failure, Door ajar

Refrigeration Type: Forced air refrigeration system

Refrigerant: R600a, CFC Free

Condenser &amp; Evaporator: Bundy tube condenser

Defrost: Auto defrost

Door: Glass door with heat reflection film

Construction: Structure Unibody Design &amp; Mono-assembly foaming (Rigid polyurethane insulation material)

Internal

Material: Embossed aluminum sheet

External

Material: PCM

Shelves 4 pcs

Door: Glass door with heat reflection film

Consumption: 136W

Power Supply: AC110/220V±10%, 50/60Hz

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE,](#)