

**Product Name :**  
Pathology Lab Autoclave**Product Code :**  
MED-PTH0056**Description :**

Pathology Lab Autoclave

**Technical Specification :**

Easy operation, safe and reliable.

The sealing ring is made of silicon rubber.

Two safety valves keep double safety (for timing control type).

It has self-adjustment function (for timing control type)

It is well sealed and effectively prevents steam leakage.

Over temperature & over pressure auto-protection.

Electric heating or gas heating can be used (for Dual fuel Type).

The sterilizer is provided with a power switch, user can control the operation easily.

Immersion type electric heating pipe, thermal efficiency is high, stainless steel material tube is durable to use.

Steam-water inner circulation system: No steam discharge and the environment for sterilizing will be clean and dry.

Double groove pressure gauge, could clearly indicate, convenient observation, reliable performance.

When at lower water level, or at the end of heat preservation time, have buzzer alarm, timing control type will automatically stop heating (for timing control type).

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**