

**Product Name :**  
Pathology Lab Mixers

**Product Code :**  
MED-PTH0058



**Description :**

Pathology Lab Mixers

**Technical Specification :**

Shaking Mode: Orbital  
Orbital Diameter: 4mm  
Motor Input: 58W  
Motor Output: 10W  
Speed Range: 0~2500rpm  
Speed Display: Scale  
Operation Mode: Continuous/touch operation  
Ambient Temp: 5-40°C  
Relative Humidity: 80%  
Protection Type: IP21  
Power Consumption: 60W  
Power supply: AC110/220V±10%, 50/60Hz

