

Product Name :
Pathology Lab Shakers**Product Code :**
MED-PTH0059**Description :**

Pathology Lab Shakers

Technical Specification :

Shaking Motion: Orbital /Linear
Orbital Diameter: 10mm
Max: Shaking Weight
(with attachment): 7.5kg
Motor Type: Brushless DCmotor
Motor Input: 28W
Motor Output: 15W
Speed Range (rpm): 100~500 (Orbital) 100~350 (Linear)
Speed Display: LCD
Time Range: 1~1199min
Operation Mode: Continuous/timed operation
Ambient Temp: 5~40°C
Protection Type: IP21
Connection Interface: RS232
Power Consumption: 30W
Power Supply: AC110/220V±10%, 50/60Hz

Medical Hospital Equipments Manufacturers

Website: [Medical Hospital Equipments Manufacturers](#),

NAUGRA[®]