

**Product Name :**  
Pathology Lab Shakers

**Product Code :**  
MED-PTH0059



**Description :**

Pathology Lab Shakers

**Technical Specification :**

Shaking Motion: Orbital /Linear  
Orbital Diameter: 10mm  
Max: Shaking Weight  
(with attachment): 7.5kg  
Motor Type: Brushless DCmotor  
Motor Input: 28W  
Motor Output: 15W  
Speed Range (rpm): 100~500 (Orbital) 100~350 (Linear)  
Speed Display: LCD  
Time Range: 1~1199min  
Operation Mode: Continuous/timed operation  
Ambient Temp: 5~40°C  
Protection Type: IP21  
Connection Interface: RS232  
Power Consumption: 30W  
Power Supply: AC110/220V±10%, 50/60Hz

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**