

Product Name :
Pathology Lab Stirrers

Product Code :
MED-PTH0060



Description :

Pathology Lab Stirrers

Technical Specification :

Work Plate Size: 120*120mm
Work Plate Material: Aluminum
Max. Stirring Capacity: 2L
Speed Range: 100-2000rpm
Heating Temp: 380 °C
Power Consumption: 180W
Power Supply: AC110/220V±10%, 50/60Hz
Standard Accessory: Adjustable rack, stirrer bar; PT100 sensor (only for BS-4HC)