

Product Name :
Pathology Lab Automatic Slide/Tissue Stainer

Product Code :
MED-PTH0066



Description :

Pathology Lab Automatic Slide/Tissue Stainer

Technical Specification :

Cup Number: 18 (Station 1 is for drying, Station 10 is for washing)
Single Cup Capacity: 750ml
Each Cup Processing Time: Adjustable within 1 hour
Slides Processed Simultaneously: 52 slides
Dripping Time: Adjustable within 0~30 s (shakes above cup)
Consumption: 500W
Power Supply: AC220V±10%, 50/60Hz; AC110V±10% 60Hz