

**Product Name :**  
Pathology Lab Automatic Slide/Tissue Stainer

**Product Code :**  
MED-PTH0066



**Description :**

Pathology Lab Automatic Slide/Tissue Stainer

**Technical Specification :**

Cup Number: 18 (Station 1 is for drying, Station 10 is for washing)  
Single Cup Capacity: 750ml  
Each Cup Processing Time: Adjustable within 1 hour  
Slides Processed Simultaneously: 52 slides  
Dripping Time: Adjustable within 0~30 s (shakes above cup)  
Consumption: 500W  
Power Supply: AC220V±10%, 50/60Hz; AC110V±10% 60Hz