

Product Name :
Pathology Lab Tissue Embedding Center**Product Code :**
MED-PTH0068**Description :**

Pathology Lab Tissue Embedding Center

Technical Specification :

Paraffin Chamber Capacity: 3L
Temp.Range of: Forceps Wells 55~70°C
Paraffin-melting Chamber
Thermal Storage Compartments
Heated Working Areas
Paraffin Dispenser
Temp. Precision: $\pm 1\%$
Material of Work Table: Granite, easy to clean
Paraffin Flow Control: Finger touch plate and optional foot pedal
Power Supply: AC110/220V $\pm 10\%$, 50/60Hz
Consumption: 650W

