

Product Name :
Low Temperature Freezer (Vertical Model)**Product Code :**
NMCBS11054**Description :**

Our company offers low temperature freezers that are designed for freezing and storing of serum, vaccines, biological and medical specimens, clinical samples etc. at low temperature. These freezers are known for their robust construction, high efficiency and long service life. Some of the features of these freezers are:

- The System includes micro controller based PID temperature controller with digital temperature display of set value and process value.
- Temperature alarm system with adjustable high - low alarm limits includes audio visual alarm on temperature deviation.
- High density CFC Free Poly Urethane Foam insulated cabinet with corrosion resistant stainless steel interior; Mild Steel powder coated exterior finish (STD Model). OR Exterior & Interior made of stainless steel 304 (GMP Model).
- Keyed door lock for improved security with magnetic door gasket.
- Full length inner glass door for clear inner view
- Open stainless steel utility baskets, Low Temperature Freezer (Vertical Model)
- Air-cooled, hermetically sealed compressor unit coupled with evaporation coil and condenser.
- Temperature range Ambient to -20°C, -40°C Low Temperature Freezer (Vertical Model), Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers
- Resolution : +/-1°C

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,