

**Product Name :**  
Industrial Refrigerators**Product Code :**  
NMCBS11056**Description :**

Industrial Refrigerators We offer a wide range of Refrigerators that provides an ideal freezing environment for advanced medical and industrial applications. These are basically used in prevention of plasma, related blood components, vaccines, micro-organism, and testing of materials and electronic components at -80 Deg C. Some of the features are:

- Construction - Double Walled Construction, Exterior made of mild steel powder coated & interior made of stainless steel 304 (STD Model) or Exterior & interior made of stainless steel 304 (GMP Model)
- High density CFC Free Poly Urethane Foam insulated cabinet
- Keyed door lock for improved security with magnetic door gasket.
- Full Length inner glass door for clear inner view.
- Environment friendly Ultra Bio Freezers with micro control and CFC free cascade refrigeration system.
- Air cooled, hermetically sealed compressor unit.
- Temperature Range: -80°C controlled by microprocessor based temperature controllers.
- Temperature alarm system with adjustable high - low alarm limits includes audio visual alarm on temperature deviation.
- 230V Ac 50 HZ Mains. Further, the technical specifications are as follows: Order Code Capacity Cu.ft Ltrs Inner Dimensions H x D x W (cms) No. of Shelves

Order Code	Capacity Cu.ft	Ltrs	Inner Dimensions H x D x W (cms)	No. of Shelves
GS - 5a	3 cft	100L	50 x 45 x 45	2
GS - 5b	6 cft	180L	70 x 51 x 51	2
GS - 5c	10 cft	300L	100 x 55 x 55	3

Industrial Refrigerators, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

