

**Product Name :**  
Ultra-Low Degree Deep Freezer Laboratory

**Product Code :**  
NMCBS11058



**Description :**

Ultra-Low Degree Deep Freezer Laboratory - Vertical Deep Freezer for **Hospital Refrigeration System**  
Exporters, **Medical Lab Refrigeration Equipments** Manufacturers

**Technical Specification :**

Ultra-Low Degree Deep Freezer Laboratory - **Vertical Deep Freezer Capacity 158L**  
Condenser & Evaporator Made of copper

**Construction Structure:**

Unibody Design & Mono-assembly foaming (Rigid polyurethane insulation material)

Internal SS304

External Material Cold-rolled steel with anti-bacteria powder coating

Door: lockable door

Caster: 4 casters

