

**Product Name :**  
Hot Air Oven Memmert**Product Code :**  
NMCBS11060**Description :**

Hot Air Oven Memmert We offer hot air laboratory ovens that are used in a variety of thermal processing applications including general lab work, component and stability testing, core hardening, drying glassware, and sterilizing. These are available in various specifications thus meeting the varied needs of our clients. Some of the features are as follows: •Double Walled – Outer MS – Inner MS or St.Steel. •80/20 Ni Heating Elements placed at two vertical sides in addition to the bottom one for uniform heating in the chamber. •Fan arrangement for forced air circulation through out the chamber. •Temp. Range : 5°C above ambient to 300°C 2°C •230V AC 50HZ Mains Supply. •Digital Temperature Controller. Models Chamber H x W x D (cms) Watts Shelves AEI-20a 30 x 30 x 30 750 1 AEI-20b 35 x 35 x 35 1000 2 AEI-20c 45 x 45 x 45 1500 2 AEI-20d 45 x 45 x 60 1750 2 AEI-20e 60 x 60 x 60 2000 3 AEI-20f 45 x 60 x 90 2200 3 AEI-20g 60 x 60 x 90 2500 3