

Product Name :
Water Baths Thermostatic**Product Code :**
NMCBS11063**Description :**

Water Baths Thermostatic We are offering a wide range of water baths thermostatic that are fabricated from top quality materials as per the international quality standards. These water baths are known for their high efficiency, robust construction and corrosion resistant. Some of the features are: •Double walled construction •Inner made of stainless steel 304 & Outer made of mild steel with powder coating finish. •Temperature controlled by Thermostat. •Temperature Range : 5°C above ambient to 95 °C •6 or 12 holes 7.5cm dia with set of concentric rings for each hole. Water Baths Thermostatic, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,