

Product Name :
Viscosity Bath**Product Code :**
NMCBS11064**Description :**

Viscosity Bath We offer a wide range of viscosity bath that are widely used for tests in the field of medical, agricultural, bio- chemical, research & development laboratories. These baths are of double walled construction, exterior made of mild steel powder coated, interior made of stainless steel 304 with puff insulation, S.S lid on top with the provision of viscometer holders. Full size toughened glass windows at front & back side of the bath with illumination arrangement. Some of the technical specifications of these water baths are: Circulation: Water circulation by AC/DC geared motor stirrer mounted on a base and fitted on the top of the bath. Heaters: S.S tubular heaters for maintain temperature. Refrigeration System: Hermetically sealed compressor provides temperature for below ambient conditions. Refrigerant: CFC free. Temperature Controller: Digital Temperature controller with PT-100 sensor. Temperature Range: 5°C above ambient to 95 °C Accuracy: ± 0.5 °C Viscosity Bath, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers