

**Product Name :**  
Autoclaves Vertical**Product Code :**  
NMCBS11068**Description :**

Autoclaves Vertical We are offering a wide assortment of autoclaves in vertical sizes. These autoclaves are manufactured from quality materials thus providing maximum satisfaction to our clients. Some of the features of these autoclaves are: •Double Walled inside made of thick stainless steel and outer body of heavy gauge mild steel attractively painted (STD Model) or interior and exterior made of stainless steel 304 (GMP Model). •Lid is of SS die pressed and fitted to the body / container tank by means of wing nuts & bolts working for fly nut system. •The lid is fitted with a pressure gauge, safety valve, steam release cock and extra safety valve. •Max pressure 20 psi / Working Pressure 15 psi. •Tested up to 40 psi. •Foot Pedal Lifting Arrangement (Optional) •Water Level Arrangement (Optional) •Perforated SS Basket (Optional) •Automatic Pressure Switch to cut off supply at set pressure (Optional). •Radial Locking System (Chrome Plated) with paddle lifting device In place of Wing nut System (Optional) Further, the technical specifications of these autoclaves are: Size Dia x Ht. Cap (Litre) Rating (KW) a) 25 x 45 22 1.5 b) 30 x 50 35 2 c) 35 x 55 53 3 d) 45 x 60 95 4 e) 45 x 75 119 4.5 f) 55 x 75 178 6 g) 55 x 90 210 8 Autoclaves Vertical, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers