

Product Name :
Low Temperature Cryostat Baths/Recirculators**Product Code :**
NMCBS11069**Description :**

Low Temperature Cryostat Baths/Recirculators Our company also offers low temperature cryostat baths that are of mild steel exterior finished with powder coating paint. These baths are made of S.S and have immersion type heaters for accuracy of temperature, with PUF insulation. We offer two types of baths i.e., Table Top and Floor Models. Model Temp. Range Capacity (Ltr.) Chamber Dimension W x D x L (cms) Bath Opening (cms) Table top models AEI-34a 0° to 50°C 12 Ltr 20 x 15 x 45 20 x 20 AEI-34b -20° to 50°C 12 Ltr 20 x 15 x 45 20 x 20 AEI-34c -40° to 50°C 12 Ltr 20 x 15 x 45 20 x 20 Floor models AEI-34d 0° to 50°C 15 Ltr 20 x 17 x 45 20 x 20 AEI-34e -20° to 50°C 15 Ltr 20 x 17 x 45 20 x 20 AEI-34f -40° to 50°C 15 Ltr 20 x 17 x 45 20 x 20 AEI-34g 0° to 50°C 25 Ltr 25 x 20 x 52 25 x 25 AEI-34h -20° to 50°C 25 Ltr 25 x 20 x 52 25 x 25 AEI-34i -40° to 50°C 25 Ltr 25 x 20 x 52 25 x 25 Low Temperature Cryostat Baths/Recirculators, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,