

Product Name :
Orbital Shaking Incubators**Product Code :**
NMCBS11070**Description :**

Orbital Shaking Incubators We are leading manufactures and suppliers of Orbital Shaking Incubators that are made up from top graded raw material. The RPM range varies between 50 to 250 RPM complete with steel tray with lotus shaped spring action clamps. Our products stand on international quality standard and are appreciated for their high efficiency, ease operation and low maintenance. Some of the features of our various models are as follows: Models Capacity (Flasks) AEI-8a 9 x 250 ml AEI-8b 16 x 250 ml AEI-8c 20 x 250 ml AEI-8d 36 x 250 ml AEI-8e 9 x 500 ml AEI-8f 16 x 500 ml AEI-8g 20 x 500 ml AEI-8h 36 x 500 ml Orbital Shaking Incubators, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,