

**Product Name :**  
Photo Stability Chambers**Product Code :**  
NMCBS11071**Description :**

Photo Stability Chambers Our company offers wide range of high quality Photo Stability Chambers that are designed with double walled SS mirror polished interior that includes air plenum, cabinet, supports and shelves. Photo stability chambers have UV tubes for ultra violet rays, integrated UV energy share of 200 watts and puff insulation fluorescent tubes for day light effect. Control System: Precise control of Temperature by using microprocessor based PID Controller with out volatile memory. Temperature Range: 15deg.C to 40 deg.C +/-0.2deg.C Uniformity : 0.5deg.C Some of the technical specifications are as follows: Model Capacity Inner Chamber H x D x W (cms) No. of Components Cu.Ft Ltrs AEI-10a 6 cft 180 Ltr 70 x 51 x 51 3 AEI-10b 10 cft 300 Ltr 100 x 55 x 55 4 AEI-10c 12 cft 360 Ltr 100 x 60 x 60 6 AEI-10d 16 cft 480 Ltr 114 x 65 x 65 Photo Stability Chambers, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers 8

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#),