

Product Name :
Dehumidifiers**Product Code :**
NMCBS11074**Description :**

Dehumidifiers We offer Dehumidifiers to control environment that are constructed to control positive humidity. Dehumidifiers are designed using stainless or mild steel which are very portable on refrigeration system. They are highly appreciated for dehumidifiers of room atmosphere for pharmaceuticals, food industry, packaging rooms & highly sophisticated instrument rooms where lower humidity is required. Some of the technical specifications of these are as follows: Model Coverage Area Height Capacity AEI-33a 25 Sq. meters 2.5 meters 1 T AEI-33b 35 sq. meters 2.5 meters 1.5 T AEI-33c 50 sq. meters 2.5 meters 2 T AEI-33d 65 sq. meters 2.5 meters 2.5 T Dehumidifiers, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,