

Product Name :
Rotary Shaker**Product Code :**
NMCBS11078**Description :**

Rotary Shaker We present a wide array of premium quality Rotary Shaker to our esteemed clients. The premium quality Rotary Shaker are platform style mounted on heavy duty MS angle with heavy Iron pulley cast. The shakers are driven electrically with PMDC motor presenting the rotary shaking stroke of preset speed of 200 rpm or variable pace of 50 – 200 rpm. Some of the technical specifications of these shakers are: Model Capacity To hold Flasks Platform Size L x B (cms) AEI-25a 9 x 500ml. 50 x 50 AEI-25b 12 x 500ml. 50 x 50 AEI-25c 16 x 500ml 60 x 60. AEI-25d 25 x 500ml 75 x 75 AEI-25e 36 x 500ml 80 x 80 AEI-25f 49 x 500ml 90 x 90 AEI-25g 100 x 500ml 120 x 120 Rotary Shaker, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,