

**Product Name :**  
Laboratory Incubator**Product Code :**  
NMCBS11080**Description :**

Laboratory Incubator We offer laboratory incubator that is used for incubation of bacteria and help in decomposing organic matter in waste water. This Incubator is also used in biology, botany, virology and oceanography. Some of the features of these incubators are: • Double Walled Construction • Exterior made of mild steel powder coated and interior made of stainless steel 304 or Exterior and Interior made of stainless steel 304. • 80/20 Ni heating elements. • Temperature Range: 5°C above ambient to 80°C • Digital Temperature controller, • Built in glass window in the front door is provided for inner observation • Power supply: 230V AC Mains 50 HZ. Further, the technical specifications of these incubators are: Models Chamber Watts Shelves H x W x D (cms) AEI-18a 30 x 30 x 30 250 1 AEI-18b 35 x 35 x 35 300 2 AEI-18c 45 x 45 x 45 500 2 AEI-18d 60 x 45 x 45 600 2 AEI-18e 60 x 60 x 60 750 3 AEI-18f 90 x 60 x 45 900 3 AEI-18g 90 x 60 x 60 1000 3 Laboratory Incubator, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,