

**Product Name :**  
Environmental Test Chamber**Product Code :**  
NMCBS11082**Description :**

Environmental Test Chamber Double walled construction Exterior made of stainless steel 304 Interior made of stainless steel 316 Puff insulation (75 mm) Full length inner glass door for clear inner view Fan arrangement for forced air circulation through out the chamber for uniform temperature Pid temperature controller cum indicator Cooling system : hermetically sealed compressor unit coupled with evaporation coil & condenser. Temperature range : 10<sup>o</sup>cto 60<sup>o</sup>c Resolution: +/-1<sup>o</sup>c Humidity range : 40%rh to 95%rh Resolution: +/-3%rh Environmental Test Chamber, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,