

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Laboratory Oven

Product Code :
NMCBS11084



Description :

Laboratory Oven Available in two different models- STD & GMP models -Exterior made of mild steel and interior made of stainless steel 304. Air Heating elements are placed three vertical side. Glass wool insulation for preventing heat loss Laboratory Oven, Heat And Refrigeration System Exporters, Heat And Refrigeration System Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,